

12-Channel Adjustable Volume Micro Pipettes



12 channel Micropipettes

Description

- Pipettes cover a volume range from 0.5µL to 300µL.
- Ergonomic design provides excellent operating experience.
- Large display window allows for easy volume identification.
- Easy calibration and maintenance.
- 12 channel pipettes are appropriate for 96 well plates.
- Dispensing head rotates for effortless pipetting convenience.
- Individual piston and tip cone assemblies allowing easy repair and maintenance.
- Compound material-made tip cone secures high sealing performance.
- Compatible with most universal tip brands.

Specifications

MSPIPET12SWISS 12-channel Adjustable Volume Pipettes				
Volume Range	Increment	Test Volume	(Accuracy error)	(Precision error)

%	µL	%	µL						
0.5-10µL		0.1µL		10µL		1.50	0.15	1.50	0.15
5µL		2.50	0.125	2.50	0.125				
1µL		4.00	0.04	4.00	0.04				
5-50µL		0.5µL		50µL		1.00	0.5	0.50	0.25
25µL		1.50	0.375	1.00	0.25				
5µL		3.00	0.15	2.00	0.1				
10-100µL		1µL		100µL		0.80	0.8	0.30	0.3
50µL		1.00	0.5	0.80	0.4				
10µL		3.00	0.3	2.00	0.2				
50-300µL		5µL		300µL		0.70	2.1	0.25	0.75
150µL		1.00	1.5	0.50	0.75				
50µL		1.50	0.75	0.80	0.4				