

Constant Temperature Drying Oven, 43L



MSBJPXHDO43SWISS Constant Temperature Drying Oven: Specification:

- **Model:** MSBJPXHDO43SWISS
- **Capacity:** 43 liters
- **Heating Method:** Bottom heating, natural convection
- **Temperature Range:** RT + 10°C to 250°C
- **Display:** LCD display
- **Temperature Uniformity:** ? Max temp. $\pm 3\%$
- **Temperature Precision:** 0.1°C
- **Timing Range:** 0 to 9999 minutes/hours
- **Inner Size (L × W × H):** 352 × 351 × 346 mm
- **Number of Shelves:** 2

- **Power Consumption:** 1500 W
- **Power Supply:** Standard AC 220V, 50/60 Hz

This drying oven provides reliable and uniform heating for various drying and sterilization processes in the lab.