

8-Channel Adjustable Volume Micro Pipettes



8 channel Micropipettes

Description

- Pipettes cover a volume range from 0.5?L to 1000?L.
- Ergonomic design provides excellent operating experience.
- Large display window allows for easy volume identification.
- Easy calibration and maintenance.
- 8 channel pipettes are appropriate for 96 well plates.
- Dispensing head rotates for effortless pipetting convenience.
- Individual piston and tip cone assemblies allowing easy repair and maintenance.
- Compound material-made tip cone secures high sealing performance.
- Compatible with most universal tip brands.

Specifications

MSPIPET8SWISS 8-channel Adjustable Volume Pipettes

Volume Range

%

0.5-10?L

5?L

1?L

5-50?L

25?L

5?L

10-100?L	
50?L	
10?L	
50-300?L	
150?L	
50?L	