

Heating Block / Dry Bath



Overview of the **MSNEC211SWISS Heating Block / Dry Bath**:

MSNEC211SWISS Heating Block / Dry Bath

Features:

- **Instant Temperature Display:** Quick view of temperature status.
- **Built-in Safety Device:** Protection from over-temperature.
- **Intelligent Algorithm:** Efficient heating control.
- **Automated Power Recovery:** Automatically resumes function after a power cut.

- **Multiple Module Options:** Allows various configurations for different applications.
- **LED Display:** Clear and precise digital display.
- **Failure Detection & Alarm:** Notifies user of errors or malfunctions with a buzzing alarm.

Technical Specifications:

- **Model:** MSNEC211SWISS
- **Temperature Range:** Room temperature + 5° to 100°
- **Temperature Stability:** $\pm 0.5^\circ$
- **Display Precision:** 0.1°
- **Heating Up Time:** ? 12 minutes (from 25° to 100°)
- **Timing Range:** 5 seconds to 23 hours 59 minutes
- **Power:** 96W
- **Power Supply:** AC100-240V/1.5A to DC24V/4A
- **Module Dimensions:** 111 x 76 x 33 mm
- **Product Size:** 171 x 161 x 103 mm

Adaptive Module Options:

- **MSEC8:** Centrifuge tube module, accommodating six 50 ml tubes
- **Custom Modules:** Available per customer requirements