

Trinocular Head Microscope



Trinocular Head Microscope

Specification

- **Model:** MSXS208BSWISS
- **Head:** Sliding Trinocular Head Inclined at 45°
- **Nosepiece:** Ball Bearing Quadruple Nosepiece
- **Eyepiece:**
 - Wide Field Plane-scope Eyepiece: WF10X/18mm
- **Objective:**
 - Achromatic 4X, 10X, 40X (Spring), 100X (Spring, Oil Immersion)
- **Stage:**
 - Double Layer Mechanical Stage (135 × 140 mm)
- **Condenser:**
 - ABBE NA 1.25 Condenser with Iris Diaphragm & Filter, Rack & Pinion

Adjustable Condenser

- **Illumination:**
 - Built-in LED Lamp
- **Power Supply:**
 - AC 110/220V $\pm 10\%$, 50/60Hz
- **Packaging:**
 - Carton with Foam
- **Package Size:** 250 × 360 × 460 mm
- **Gross Weight:** 8 kg