

Vortex Mixer



Key Features:

1. **Compact Design:** Small, lightweight, and easy to carry for convenient usage.
2. **Brushless DC Motor:** Provides strong, maintenance-free operation for reliable performance.
3. **Durable Silicone Vibrating Head:** Imported silicone head with high toughness, flexibility, and long service life.
4. **Fast Mixing:** Can swirl samples within 5 seconds for quick processing.
5. **Press-to-Start with Stepless Speed Control:** Easy operation with adjustable speed settings, suitable for various applications.
6. **Stable Operation:** Features flexible anti-vibration foot pads, a metal base, and a triangular structure for enhanced stability, low noise, and no displacement during operation.

Product Specifications:

- **Model:** MSLCHMXM2SWISS

- **Voltage:** DC 12V 1A
- **Power:** 5W
- **Speed Range:** 500-3500 rpm
- **Drive Mode:** Brushless DC motor
- **Operation Mode:** Press and jog operation
- **Circumference Diameter:** 4.5 mm
- **Maximum Processing Capacity:** 50 ml
- **Tube Diameter:** Max 35 mm
- **Net Weight:** 540 g
- **Dimensions:** 110 × 110 × 75 mm
- **Gross Weight:** 780 g
- **Package Dimensions:** 210 × 155 × 135 mm

This Vortex Mixer is designed for quick, efficient sample mixing with durable construction and easy-to-use features, making it ideal for laboratory environments.