

Water Distiller, Water Control Type



Distilled Water Machine

This machine is ideal for producing distilled water in laboratories within the medical and health industry, chemical industry, and scientific research units.

High-Quality Construction: Made from premium stainless steel using advanced manufacturing techniques.

Safety Features: Equipped with water-deficient and power-off alarms.

Automatic Functionality: Automatically replenishes water and reheats.

Operational Indicators: Includes power heating instructions and water level display.

Electrode-Free Principle: Prevents power-off issues due to electrode fouling from long-term use.

Specifications

Model: MSWDST5ZLSWISS

Water Output: 5L/h

Power: 5kW

Electricity: AC220V 50Hz

Overall Dimension (D*H): 310mm × 740mm

Net Weight/Gross Weight: 7.1kg / 7.8kg