

## 12-Channel Adjustable Volume Micro Pipettes



### 12 channel Micropipettes

#### Description

- Pipettes cover a volume range from 0.5µL to 300µL.
- Ergonomic design provides excellent operating experience.
- Large display window allows for easy volume identification.
- Easy calibration and maintenance.
- 12 channel pipettes are appropriate for 96 well plates.
- Dispensing head rotates for effortless pipetting convenience.
- Individual piston and tip cone assemblies allowing easy repair and maintenance.
- Compound material-made tip cone secures high sealing performance.
- Compatible with most universal tip brands.

#### Specifications

<b>MSPINET12SWISS 12-channel Adjustable Volume Pipettes</b>				
Volume Range	Increment	Test Volume	(Accuracy error)	(Precision error)

%	?L	%	?L						
0.5-10?L		0.1?L		10?L		1.50	0.15	1.50	0.15
5?L	2.50	0.125		2.50		0.125			
1?L	4.00	0.04		4.00		0.04			
5-50?L		0.5?L		50?L		1.00	0.5	0.50	0.25
25?L	1.50	0.375		1.00		0.25			
5?L	3.00	0.15		2.00		0.1			
10-100?L		1?L		100?L		0.80	0.8	0.30	0.3
50?L	1.00	0.5		0.80		0.4			
10?L	3.00	0.3		2.00		0.2			
50-300?L		5?L		300?L		0.70	2.1	0.25	0.75
150?L	1.00	1.5		0.50		0.75			
50?L	1.50	0.75		0.80		0.4			