


-25? Freezer, 100L**-25? Freezer, 100L****Specification**

- **Model:** MSBDF25V100SWISS
- **Capacity:** 100L
- **Type:** Vertical
- **Temperature Range:** -10? to -25?
- **Temperature Accuracy:** 0.1?
- **Control System:** Microprocessor control with a digital display
- **Alarm System:**
 - High and low temperature alarm
 - Sensor failure alarm
 - Door open alarm
- **Refrigeration Type:** Direct refrigeration

- **Refrigerant:** CFC-free
- **Condenser:** Bundy tube