


## Dry Bath



### Model: MSDBI1001SWISS

- **Standard:** 2 blocks
- **Temperature Range:** RT.+5°C ~ 160°C
- **Temperature Accuracy:** ±1°C (@120°C)
- **Temperature Uniformity:** ±0.5°C
- **Display Accuracy:** 0.1°C
- **Time Range:** 1-99h59min / ?
- **Heating Time:** ?15 min (from 25°C to 160°C)
- **Max Power:** 500W
- **Voltage:** AC 220V, 50Hz
- **Size:** 260 × 220 × 100 mm
- **Package Size:** 360 × 310 × 230 mm
- **Gross Weight:** 6 kg

