


Electrophoresis Power Supply



MSBEP300ISWISS Electrophoresis Power Supply: Specifications:

- **Model:** MSBEP300ISWISS
- **Voltage Output Range:**
 - 2-300V (Resolution: 1V)
- **Current Output Range:**
 - 5-2000mA (Resolution: 1mA)
- **Power Output Range:**
 - 1-300W (Resolution: 1W)
- **Input Voltage:**
 - 220V, 50Hz
- **Parallel Output:**
 - 4 groups
- **Timing Function:**

- 1 minute to 99 hours 59 minutes, adjustable

This power supply is designed to offer precise control for electrophoresis applications with a wide voltage, current, and power output range, making it suitable for various laboratory needs. The inclusion of parallel outputs for up to 4 groups provides enhanced flexibility for multiple runs.