

Large Scale Electronic Platform Balance, 150Kg



MSBE150000LSWISS Large Scale Electronic Platform Balance Specification:

- **Model:** MSBE150000LSWISS
- **Capacity:** 0 to 150,000 g
- **Readability:** 1 g
- **Scale Size:** 500 × 400 mm
- **Power Supply:** 100~240V, 50/60Hz
- **External Size (W × D × H):** 530 × 300 × 680 mm
- **Net Weight:** 11 kg
- **Package Size (W × D × H):** 870 × 465 × 225 mm
- **Gross Weight:** 14 kg

