


Laboratory Incubator, 30L



MSBJPXH30LSWISS Laboratory Incubator: Specification:

- **Model:** MSBJPXH30LSWISS
- **Capacity:** 30 liters
- **Temperature Control Method:** PT100
- **Temperature Range:** RT + 5°C to 85°C
- **Temperature Fluctuation:** ±0.2°C (at 37°C)
- **Temperature Uniformity:** ±0.5°C (at 37°C)
- **Circulation Mode:** 4 sides heating, forced convection
- **Ambient Temperature:** 5°C to 30°C (recommended: 25°C ± 2°C)
- **Timing Range:** 0 to 9999 minutes/hours
- **Shelves:** 2 pcs

- **Power Consumption:** 150 W
- **Display:** LCD
- **Power Supply:** AC 220V \pm 10%, 50/60 Hz (Standard)
- **Internal Size (L x W x H):** 312 x 314 x 306 mm

This incubator is designed for precise temperature control and uniform heating, suitable for various laboratory applications.