

## Trinocular Head Microscope



### Trinocular Head Microscope

#### Specification

- **Model:** MSXS208BSWISS
- **Head:** Sliding Trinocular Head Inclined at 45°
- **Nosepiece:** Ball Bearing Quadruple Nosepiece
- **Eyepiece:**
  - Wide Field Plane-scope Eyepiece: WF10X/18mm
- **Objective:**
  - Achromatic 4X, 10X, 40X (Spring), 100X (Spring, Oil Immersion)
- **Stage:**
  - Double Layer Mechanical Stage (135 × 140 mm)
- **Condenser:**
  - ABBE NA 1.25 Condenser with Iris Diaphragm & Filter, Rack & Pinion

Adjustable Condenser

- **Illumination:**
  - Built-in LED Lamp
- **Power Supply:**
  - AC 110/220V  $\pm 10\%$ , 50/60Hz
- **Packaging:**
  - Carton with Foam
- **Package Size:** 250 × 360 × 460 mm
- **Gross Weight:** 8 kg