

Moisture Analyzer



Moisture Analyzer

MSMBA1102SWISS Halogen Moisture Analyzer

Key Features:

- 1. High-Brightness LED Display:** Equipped with a clear and easy-to-read LED display screen, ensuring visibility even in various lighting conditions.
- 2. Durable Aluminum Sample Weighing Pan:** The pan is solid, durable, and easy to clean, providing long-term reliability and ease of maintenance.
- 3. Double Bond Design:** Facilitates easy operation, making the analyzer user-friendly and efficient.
- 4. Comprehensive Protection Measures:** Includes overload protection, over-temperature protection, door lock protection, and multiple other protective features to ensure safe and reliable operation.
- 5. Built-In Temperature and Time Settings:** Suitable for various types of sample testing with customizable parameters, offering one-key operation for convenience.
- 6. High-Performance Halogen Heating Lamp:** Provides stable and quick heating, ensuring efficient and consistent moisture analysis.

Technical Specifications:

- **Model:** MSMBA1102SWISS
- **Capacity:** 110g
- **Readability:** 0.01g
- **Tray Size:** 96mm
- **Heating Mode:** Halogen Lamp
- **Display:** 5-inch touch screen
- **Temperature Setting Range:** 40~160
- **Measuring Range:** 90%~10%
- **Storage Capacity:** 20 sets of data
- **Communication Interface:** RS232
- **Printer:** Optional

This summary highlights the key features and technical specifications of the MSMBA1102SWISS Halogen Moisture Analyzer, emphasizing its reliability, user-friendly design, and suitability for precise moisture analysis in laboratory settings.