


Static Storage Liquid Nitrogen Container



Here are the details for the **Static Storage Liquid Nitrogen Container (Model: MSYDS35S6SWISS)**:

Technical Specifications:

- **Capacity:** 35.5 liters
- **Caliber (Opening Diameter):** 50 mm
- **Accessories:**
 - **Cannisters:** 6 units for sample storage
 - **Lockable Lid:** 1 unit for secure closure
 - **Protective Sleeve:** 1 unit for added protection during storage and handling

Key Features:

- 1. High Capacity:** With a 35.5-liter storage capacity, it is suitable for long-term storage of biological samples.
- 2. Compact Opening:** The 50 mm caliber allows for effective insulation and preservation of liquid nitrogen.
- 3. Cannister Design:** Six cannisters provide organized storage for biological samples.
- 4. Lockable Lid:** Ensures security and safety of the stored contents.
- 5. Protective Sleeve:** Designed to enhance durability and safety during transport and storage.

This liquid nitrogen container is ideal for laboratories or medical facilities that require secure, long-term static storage of samples at cryogenic temperatures.