

## Ultrasonic Bath /Ultrasonic Cleaner, 10L



**Premium Digital Ultrasonic Cleaner – High-Precision Cleaning with Advanced Features**

### **Key Features & Benefits**

? **Ultra-Low Noise & High Efficiency** – Advanced sound insulation technology ensures quiet operation, ideal for home and professional use.

? **Precision Cleaning with Adjustable Power & Time** – Set ultrasonic power from **10% to**

**100%** (in 10% increments) and customize cleaning cycles from **1 to 99 minutes** for optimal results.

? **Large LCD Display for Easy Monitoring** – Clearly view **time, temperature, and power level** for precise control.

? **Intelligent Sweep & Degas Functions** – Improves cleaning efficiency by removing gas bubbles and optimizing cavitation performance.

? **Energy-Saving Sleep Mode & Memory Function** – Automatically saves your last-used settings and enters sleep mode when inactive to conserve energy.

? **Durable Stainless Steel Construction** – Built to last with a **stainless steel tank, body, cover, and cleaning basket** for enhanced durability and corrosion resistance.

? **Integrated Heating Function** – Built-in heater reaches up to **60°C** (optional **80°C** available on request) to enhance cleaning effectiveness.

? **User-Friendly Design** – Comes with a **convenient drain valve** for easy water removal and a **sturdy carrying handle** for effortless transport.

## Technical Specifications

Model	MSUSC1040MSWISS
Max Ultrasonic Power (W)	240
Power Adjustment	10% - 100% (in 10% steps)
Max Heating Power (W)	500
Tank Capacity (L)	10
Tank Dimensions (L x W x H) mm	300 x 240 x 150
External Dimensions (L x W x H) mm	330 x 270 x 310
Number of Transducers	4

Model	MSUSC1040MSWISS
Ultrasonic Frequency (kHz)	40
Max Temperature Setting (°C)	Default: 60°C (Optional: 80°C – free upon request)
Timer Range (Minutes)	1 – 99
Cleaning Basket	Yes
Drain Valve	Yes
Carrying Handle	Yes

### Why Choose Our Ultrasonic Cleaner?

Our **Digital Ultrasonic Cleaner** is the perfect solution for deep and precise cleaning of delicate items, including jewelry, medical instruments, automotive parts, and more. With advanced **adjustable settings, high-powered transducers, and durable stainless steel construction**, this unit ensures **powerful cleaning with minimal noise**. The integrated heating function further enhances cleaning performance, making it a must-have for both professional and home use.

? **Order now** and experience **superior ultrasonic cleaning** with advanced features and unmatched efficiency!