

Ultrasonic Bath /Ultrasonic Cleaner, 30L



Ultrasonic Water Bath - Model: MSVTD4151ESWISS

Description:

- **Superior Cleaning Performance:** Equipped with a high-power transducer, this ultrasonic water bath ensures exceptional cleaning efficiency with visible results.
- **Digital LED Display:** Provides a precise display of time and temperature settings for controlled operations.
- **Adjustable Time and Heating:** Allows flexible time settings from 1 to 30 minutes and a maximum heating temperature of 80°C, enabling a variety of cleaning processes.
- **Reliable and Sustainable:** Designed for continuous operation with super water performance for enhanced safety.
- **Durable Construction:** Crafted from high-quality stainless steel for the body, cover, and base, ensuring extended durability.

- **Versatile Cleaning Fluids:** Compatible with tap water, alcohol, organic solvents, and other fluids, making it adaptable to various cleaning applications.

Specifications:

- **Model:** MSVTD4151ESWISS
- **Ultrasonic Power:** Max 600 W
- **Heating Power:** Max 500 W
- **Tank Capacity:** Max 30 L
- **Tank Dimensions (L x W x H):** 501 x 300 x 200 mm
- **External Dimensions (L x W x H):** 531 x 330 x 390 mm
- **Transducers:** 10 pieces
- **Ultrasonic Frequency:** 40 kHz
- **Temperature Range:** Up to 80°C
- **Timer Range:** 1 to 30 minutes
- **Additional Features:**
 - Drain Valve: Yes
 - Carrying Handle: Yes
 - Cleaning Basket: Yes

Applications: This ultrasonic water bath is ideal for industries requiring high-performance cleaning of sensitive or intricate items, providing effective and safe results across various applications.