

Vertical Autoclaves, 120L



Vertical Autoclave – Model MSCLRC231JJSSWISS

Application:

This vertical autoclave is ideal for sterilizing laboratory equipment, medical waste, and various materials in food processing and inspection sectors. It uses gravity displacement to efficiently expel cold air from the sterilization chamber, ensuring effective sterilization.

Features:

- Durable stainless steel construction for long service life.
- Equipped with a transparent viewing window to monitor the sterilization process.
- A pressure gauge is provided to monitor the internal pressure during the cycle.
- A digital display shows the current temperature and time.
- A digital display shows the current temperature and time.
- A digital display shows the current temperature and time.

Technical Specifications:

- Model MSCLRC231JJSSWISS
- Capacity: 120L
- Temperature range: 121°C to 134°C
- Pressure range: 1.05 bar to 1.2 bar
- Power: 2000W
- Voltage: 230V/50Hz

Standard and Optional Accessories:

- **Basket, mesh bottom basket, tray, bucket**

This vertical autoclave offers reliable sterilization with advanced monitoring features, ideal for laboratories and medical facilities where precision and efficiency are key.