

## Water Distiller, Water Control Type



### Distilled Water Machine

This machine is ideal for producing distilled water in laboratories within the medical and health industry, chemical industry, and scientific research units.

**High-Quality Construction:** Made from premium stainless steel using advanced manufacturing techniques.

**Safety Features:** Equipped with water-deficient and power-off alarms.

**Automatic Functionality:** Automatically replenishes water and reheats.

**Operational Indicators:** Includes power heating instructions and water level display.

**Electrode-Free Principle:** Prevents power-off issues due to electrode fouling from long-term use.

---

## **Specifications**

Model: MSWDST5ZLSWISS

Water Output: 5L/h

Power: 5kW

Electricity: AC220V 50Hz

Overall Dimension (D\*H): 310mm × 740mm

Net Weight/Gross Weight: 7.1kg / 7.8kg