Ultrasonic Bath /Ultrasonic Cleaner, 30L

Model **MSVTD4151ESWISS** | Brand **Master Science** | Unit **Pc**



This ultrasonic water bath is ideal for industries requiring high-performance cleaning of sensitive or intricate items, providing effective and safe results across various applications.

KEY HIGHLIGHTS

- Versatile Cleaning Fluids

Request a Quotation
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

Key Features

Crafted from high-quality stainless steel for the body, cover, and base, ensuring extended durability.

Versatile Cleaning Fluids
Compatible with tap water, alcohol, organic solvents, and other fluids, making it adaptable to various cleaning applications.

Technical Specifications

PARAMETER	SPECIFICATION
Ultrasonic Water Bath - ModelMSVTD4151ESWISS	
Ultrasonic PowerMax 600 W	
Heating PowerMax 500 W	
Tank CapacityMax 30 L	
Tank Dimensions (L x W x H)501 x 300 x 200 mm	
External Dimensions (L x W x H)531 x 330 x 390 mm	
Transducers10 pieces	
Ultrasonic Frequency40 kHz	
Temperature RangeUp to 80°C	
Timer Range1 to 30 minutes	
Superior Cleaning PerformanceEquipped with a high-power transducer, this ultrasonic water bath ensures exceptional cleaning efficiency with visible results.	
Digital LED DisplayProvides a precise display of time and temperature settings for controlled operations.	

PARAMETER	SPECIFICATION
Adjustable Time and HeatingAllows flexible time settings from 1 to 30 minutes and a maximum heating temperature of 80°C, enabling a variety of cleaning processes.	
Reliable and SustainableDesigned for continuous operation with super water performance for enhanced safety.	
Drain Valve	Yes
Carrying Handle	Yes
Cleaning Basket	Yes