# Water Distiller, Water Control Type

Brand **Master Science** | Unit **Pc**



This machine is ideal for producing distilled water in laboratories within the medical and health industry, chemical industry, and scientific research units.

KEY HIGHLIGHTS

- Automatic Functionality
- Operational Indicators
- Electrode-Free Principle

**Request a Quotation**
Contact us for pricing, configuration options and delivery details.

support@themsci.com
www.themsci.com

Scan to view online

## Key Features

EQUIPPED WITH WATER-DEFICIENT AND POWER-OFF ALARMS.

**Automatic Functionality**
Automatically replenishes water and reheats.

**Operational Indicators**
Includes power heating instructions and water level display.

**Electrode-Free Principle**
Prevents power-off issues due to electrode fouling from long-term use.

## Technical Specifications

| PARAMETER | SPECIFICATION |
|---|---|
| Water Output | 5L/h |
| Power | 5kW |
| Electricity | AC220V 50Hz |
| Overall Dimension (D*H) | 310mm × 740mm |
| Net Weight/Gross Weight | 7.1kg / 7.8kg |